```
第一步、判断用户权限组权限(linq) (Category = 系统 或 模块, 111 影像诊断) :
var q =
   from Customer Group Mapping c
   from Permission Group Mapping g
   from RolePermission p
   where c.CustomerId == 10 && p.SystemName == "报告书写" && p.Category
== 111 && c.GroupId == g.GroupId && g.PermissionId == p.Id
   select new
       g.Config
   };
第二步、判断是否有权限
var r = q.FirstOrDefault();
if (r == null) //没有权限
{
   return false;
}
第三步、判断是否有检查类型CT的书写权限
Config Json内容如下:
  "LimitRange": "CR,CT,MR"
var j = JsonHelper.ParseFormByJson < LimitRangeConfig > (r.Config);
if (j == null) //没有权限
{
    return false;
}
var limitRange = "," + j.LimitRange + ",";
if (limitRange.IndexOf("," + "CT" + ",", 0) < 0) //没有权限
{
   return false;
}
return true; //权限验证通过
```

# 附录 (几种常见场景Sql):

### 1、判断用户是否有相应权限

# SELECT g.config

FROM Customer\_Group\_Mapping c, Permission\_Group\_Mapping g, RolePermission p

WHERE c.CustomerId = 18 AND p.SystemName = "报告书写" AND p.Category = 111 AND c.GroupId = g.GroupId AND g.PermissionId = p.Id

### 2、获取有相应权限的所有用户

#### SELECT a.\*

FROM Customer a, Customer\_Group\_Mapping c, Permission\_Group\_Mapping g, RolePermission p

WHERE a.ld = c.Customerld AND p.SystemName = "报告书写" AND p.Category = 111 AND c.GroupId = g.GroupId AND g.PermissionId = p.Id

# 3、获取用户的所有权限

SELECT p.Id, p.RoleName, p.SystemName, p.Kind, p.Category,p.Module FROM Customer\_Group\_Mapping c, Permission\_Group\_Mapping g, RolePermission p

WHERE c.CustomerId = 18 AND c.GroupId = g.GroupId AND g.PermissionId = p.Id

#### 4、获取服务中心的所有默认用户组和权限

SELECT p.Id, p.RoleName, p.SystemName, p.Kind, p.Category, p.Module, p.DefaultGroupName, p.DefaultGroupKind FROM RolePermission p
WHERE DefaultGroupKind = 10 or DefaultGroupKind = 110 or DefaultGroupKind = 111

#### 5、获取机构的所有默认用户组和权限

SELECT p.Id, p.RoleName, p.SystemName, p.Kind, p.Category, p.Module, p.DefaultGroupName, p.DefaultGroupKind FROM RolePermission p
WHERE DefaultGroupKind = 100 or DefaultGroupKind = 110 or DefaultGroupKind = 111

## 6、获取平台的所有默认用户组和权限

SELECT p.Id, p.RoleName, p.SystemName, p.Kind, p.Category, p.Module, p.DefaultGroupName, p.DefaultGroupKind FROM RolePermission p
WHERE DefaultGroupKind = 1 or DefaultGroupKind = 11 or DefaultGroupKind = 111